

A new UV photocatalytic dosimeter based upon I_3^- formation shows that gamma irradiation and photocatalysis produce similar destructive chemical species.

Bret O'Brien, Lane D'Alessandro, Jordan Brown, Ju Chou, Patricia D. Barreto, Jose C. Barreto

Dept. of Chemistry and Mathematics, Green Technology Research Group, Florida Gulf Coast University, 10501 FGCU Blvd South, Fort Myers, FL 33965

Ionizing radiation (gamma/x-rays) produces radicals, oxidants, and reductants (ROR), useful for decontamination, which can be measured by dosimetry. Proposed dosimeters for gamma rays are based upon ROR mediated formation of tri-halide anions. We have developed a photocatalytic dosimeter that uses safer, long-wave UV irradiation, measuring tri-iodide (I_3^-) product formation from iodide (I^-). A KI solution in buffer (pH 7.2) was mixed with the commercial TiO_2 photocatalyst Evonik P-25 at a concentration of 0.1mg/mL and exposed to UV irradiation for 5 minutes. The production of I_3^- was measured by its absorbance at 350nm. The reaction kinetics were linear (plot of absorbance versus time), revealing reaction rate as the slope: 0.107 for Evonik P-25. Controls showed that if UV irradiation, the photocatalyst, or the iodide was absent, tri-iodide formation was almost not observed. We synthesized and screened several doped TiO_2 photocatalysts and 1% silver (0.067 slope) was the most effective for tri-iodide production. Our experiments demonstrate that photocatalytic oxidation of iodide can be used as a dosimeter to rank the effectiveness of photocatalysts, and that photocatalysis mimics gamma irradiation.